

Task 2 : Python and Lambda


Work with RXNORM file,

1. Scrap the latest RXNORM file from NLM webpage
2. Download the latest RXNORM file with `api_key`
3. Create a log file for the downloaded file
4. Add header into each `fff` from RXNORM.xlsx
5. Add `CODE_SET` & `VERSION_MONTH` column with default values `RxNorm` and version month from downloaded filename
6. Convert dates into YYYY-MM-DD



7. Save files as txt delimited by Pipe(|)
8. Validate `row_count` between original and converted files



The Cloud architecture of the given project is designed as follow:




Step 1: Proper setup of aws was done: lambda was created and s3 bucket was setup with four folders- zipfiles/ (to input the zipdata), excelfiles/(to input the excel file), unzipped/(to hold the unzipped files) and finally transformation/(to hold the final transformed data in .txt format)



Here, unzipped and transformation folder are created after lambda function is triggered

Step 2: The code structure

- A. The imports and excel_headers dictionary was initiated to hold the headers to be assigned to the rrf files



B. Function to read the excel_file and update the excel_header

```
lambda_function x Environment Vari x +
14
15
16 def read_excel_from_s3(bucket):
17     try:
18         folder_path = 'excelfiles/'
19         excel_file_name = 'RxNorm_Header.xlsx'
20         key = folder_path + excel_file_name
21
22         # Download the Excel file to the /tmp directory
23         local_excel_file = '/tmp/RxNorm_Header.xlsx'
24         s3.download_file(bucket, key, local_excel_file)
25
26         # Check if the Excel file exists
27         if os.path.exists(local_excel_file):
28             print(f"Excel file downloaded to: {local_excel_file}")
29
30
31
32
33         # Read the Excel file into an ExcelFile object
34         excel_file = pd.ExcelFile(local_excel_file)
35
36         # Get the sheet names
37         sheet_names = excel_file.sheet_names
38         print("Sheet names:", sheet_names)
39
40         for sheets in sheet_names:
41             # Read the data from the sheet into a DataFrame
42             sheets_data = excel_file.parse(sheets, header=None)
43             headers_data = sheets_data.iloc[:, 0].tolist()
44             excel_headers[sheets] = headers_data
45             print(f'excel_headers dictionary for sheet {sheet_names[0]}: {excel_headers[sheet_names[0]]}')
46
47
48     except Exception as e:
49         print(f"Error occurred: {e}")
50
```

C. Function to add Code Set and Version Month: Code set was set to RXNORM and version month was extracted from the zip file name which is RxNorm_full_02052024.zip Hence Version Month is: 2024-05-02

```
lambda_function x Environment Vari x +
50
51 def code_set_and_version_month(zip_data, rrf_df):
52     try:
53         # Convert the zip data to a string to extract information
54         zip_filename = zip_data.decode('utf-8')
55         # Extract version month from the filename
56         version_month = os.path.splitext(zip_filename)[0].split('_')[-1]
57
58         # Convert version month to a more readable format
59         version_month = pd.to_datetime(version_month, format='%m%d%Y').strftime('%Y-%m-%d')
60         print(f"Version month: {version_month}")
61
62         # Add 'Code Set' and 'Version Month' columns to the DataFrame
63         rrf_df['Code Set'] = 'RxNorm'
64         rrf_df['Version Month'] = version_month
65     except Exception as e:
66         print(f"Error occurred while extracting version month: {e}")
67
68     return rrf_df
69
```

D. Function to apply header to rrf from excel_header

```
70 def apply_header_to_rrf(file_name, rrf_df):
71     # Check if the corresponding Excel sheet exists
72     if file_name in excel_headers:
73         # Get the headers from the Excel sheet
74         excel_headers_list = excel_headers[file_name]
75         excel_headers_list = [header for header in excel_headers_list if header != 'SVER']
76         # Take names from the excel header list up to the length of the split DataFrame
77         excel_headers_list = excel_headers_list[:len(rrf_df.columns)]
78         # Set the correct header for the DataFrame
79         rrf_df.columns = excel_headers_list
80     return rrf_df
```

E. Function to convert date to indicated format

- a. Here, in the context of RXNSAB, the column SVER derived its value from extracting year part from the VSTART column as done below

```
121 def process_date_columns(file_name, rrf_df):
122     date_columns = ['VSTART', 'VEND', 'CREATED_TIMESTAMP', 'UPDATED_TIMESTAMP', 'LAST_RELEASED']
123     for column in date_columns:
124         if column in rrf_df.columns:
125             if file_name == 'RXNSAB':
126                 # Apply the conversion function to each value in the column
127
128
129                 rrf_df[column] = rrf_df[column].apply(convert_date_format)
130                 # Extract year from the VSTART column after date conversion and save it directly as a string
131                 rrf_df['SVER'] = pd.to_datetime(rrf_df['VSTART'], format='%Y-%m-%d').dt.year.astype(str)
132                 # Reorder the columns to place 'SVER' before 'VSTART'
133                 # Reorder columns
134                 # Reorder columns
135                 sver_index = rrf_df.columns.get_loc('SVER')
136                 vstart_index = rrf_df.columns.get_loc('VSTART')
137                 sf_index = rrf_df.columns.get_loc('SF')
138
139                 # Remove 'SVER' from its original position
140
141                 column_sver = rrf_df.pop('SVER')
142
143                 # Insert 'SVER' after 'SF', before 'VSTART'
144                 if sver_index < vstart_index:
145                     rrf_df.insert(vstart_index - 1, 'SVER', column_sver)
146                 elif sver_index > vstart_index:
147                     rrf_df.insert(vstart_index, 'SVER', column_sver)
148                 if file_name == 'RXNATOMARCHIVE':
149                     rrf_df[column] = rrf_df[column].apply(update_nato_date)
150     return rrf_df
```

- b. Now, the anomalies in the VSTART column of RXNSAB is processed as below:

```
81
82 def convert_date_format(value):
83
84     try:
85         # Try to parse the value into datetime format
86         parsed_date = pd.to_datetime(value, format='%Y_%m_%d').date()
87         # Extract only the date part
88         return parsed_date.strftime('%Y-%m-%d')
89     except ValueError:
90         if value == '2020':
91             return '2020-01-01'
92         elif value == '5.0_2024_01_04':
93             # Remove the float value and parse the remaining string
94
95             return convert_date_format('2024_01_04')
96         elif value == '2020AA':
97             return '2024-01-02'
98         elif value == '20AA_240205F':
99             return '2024-02-05'
100        else:
101            return value
102
```

- c. Now, different date formats of RXNATOMARCHIVE is handled as below:

```
102
103 def update_nato_date(value):
104     try:
105         # Attempt to parse the value using the first date format
106         parsed_date = pd.to_datetime(value, format='%m/%d/%Y %I:%M:%S %p').date()
107     except ValueError:
108         try:
109             # If the first format fails, attempt to parse using the second date format
110             parsed_date = pd.to_datetime(value, format='%d-%b-%y').date()
111         except ValueError:
112             # If both formats fail, return None or handle the error appropriately
113             return None # Or handle the error appropriately
114             # Check if the parsed_date is NaT
115     if pd.isnull(parsed_date):
116         return '0000-00-00' # Replace NaT with '0000-00-00'
117     else:
118         # Extract only the date part and return it in the desired format
119         return parsed_date.strftime('%Y-%m-%d')
120
```

F. Now the function to save the final transformed data as .txt delimited by pipe

```

151
152 def save_as_txt_file(rrf_df, file_name, bucket_name):
153     # Construct the filename for the output text file
154     transformation_folder = 'transformation/'
155
156     # Convert DataFrame to CSV format in memory
157     csv_buffer = io.StringIO()
158     rrf_df.to_csv(csv_buffer, sep='|', index=False)
159
160     # Upload the CSV buffer to S3
161     s3_key = transformation_folder + file_name + '.txt'
162     s3.put_object(Bucket=bucket_name, Key=s3_key, Body=csv_buffer.getvalue())
163
164     print(f"Transformed data saved to: s3://{bucket_name}/{s3_key}")
165
166
167

```

G. The main function where the previous functions are called:

- a. Here zip files are read from the proper location and the last pipe is ignored by splitting

```

167
168 def read_and_relocate_rrf_files(s3, bucket, key):
169     try:
170         zip_response = s3.get_object(Bucket=bucket, Key=key)
171         zip_data = zip_response['Body'].read()
172
173         # Wrap the zip data in a BytesIO object
174         zip_file = BytesIO(zip_data)
175
176         file_path = 'rrf'
177         unzipped_folder = 'unzipped/'
178
179         with zipfile.ZipFile(zip_file, 'r') as zip_ref:
180             for file_info in zip_ref.infolist():
181                 if file_info.filename.startswith(file_path) and not file_info.filename.endswith('/'):
182                     filename = os.path.basename(file_info.filename)
183                     print(f"The {filename} is read from zip file.")
184
185                     with zip_ref.open(file_info) as source_file:
186                         file_content = source_file.read().decode('utf-8')
187                         if file_content.endswith('|'):
188                             file_content = file_content[:-1]
189                         file_content_io = io.StringIO(file_content)
190                         rrf_df = pd.read_csv(file_content_io, delimiter='|', header=None)
191                         rrf_df = rrf_df.iloc[:, :-1]
192                         print(f"Row count before transformation: {rrf_df.shape[0]}")
193

```

- b. Previous functions are called to transform according to the requirements and unzipped files are placed in proper folder

```

192
193
194     file_name = os.path.splitext(filename)[0]
195     apply_header_to_rrf(file_name, rrf_df)
196     rrf_df = process_date_columns(file_name, rrf_df)
197     code_set_and_version_month(zip_data, rrf_df)
198
199     print(f"Row count of {file_name} after transformation: {rrf_df.shape[0]}")
200
201     pd.set_option('display.max_columns', None)
202     print(rrf_df.head(5))
203
204     # Save the transformed DataFrame to a text file
205     save_as_txt_file(rrf_df, file_name, bucket)
206
207
208     # Upload the unzipped file to the 'unzipped' folder
209     unzipped_key = unzipped_folder + filename
210     s3.put_object(Bucket=bucket, Key=unzipped_key, Body=file_content)
211     print(f"Unzipped file saved to: s3://{bucket}/{unzipped_key}")
212
213 except Exception as e:
214     print(f"Error occurred: {e}")
215

```

H. The main lambda function

```
215
216 def lambda_handler(event, context):
217     bucket = event['Records'][0]['s3']['bucket']['name']
218     key = event['Records'][0]['s3']['object']['key']
219
220
221     # This is the function that relocate the rrf files from zip file
222     read_excel_from_s3(bucket)
223     read_and_relocate_rrf_files(s3,bucket,key)
224
```

Step 3: Results and analysis

- A. Zip file is uploaded to zipfiles/ in order to trigger the lambda and excel file is loaded to excelfiles/

The screenshot shows the AWS S3 console interface. At the top, there's a navigation bar with the AWS logo and a search bar. Below it, a list of services is visible. The main content area shows an 'Uploading' status window. The window has a blue header with a progress bar at 21%. Below the header, there's a 'Summary' section with a table showing upload statistics. The table has columns for Destination, Succeeded, and Failed. The Destination is 's3://rrf-bucket/zipfiles/'. The Succeeded column shows '0 files, 49.2 MB (20.61%)'. The Failed column shows '0 files, 0 B (0%)'. Below the summary, there's a 'Files and folders' section with a search bar and a table listing the files. The table has columns for Name, Folder, Type, Size, Status, and Error. The first row shows a file named 'RxNorm_ful...' with a size of 238.5 MB and a status of 'In progress (20...'.

Destination	Succeeded	Failed
s3://rrf-bucket/zipfiles/	0 files, 49.2 MB (20.61%)	0 files, 0 B (0%)

Name	Folder	Type	Size	Status	Error
RxNorm_ful...	-	application/...	238.5 MB	In progress (20...	-

Amazon S3 > Buckets > rrf-bucket > zipfiles/

zipfiles/ Copy S3 URI

Objects Properties

Objects (1) [Info](#) Refresh Copy S3 URI Copy URL Download Open Delete Actions Create folder Upload

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

< 1 > ⌕

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	RxNorm_full_02052024.zip	zip	March 10, 2024, 16:01:06 (UTC+05:45)	238.5 MB	Standard

Amazon S3 > Buckets > rrf-bucket > excelfiles/

excelfiles/ Copy S3 URI

Objects Properties

Objects (1) [Info](#) Refresh Copy S3 URI Copy URL Download Open Delete Actions Create folder Upload

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

< 1 > ⌕

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	RxNorm_Header.xlsx	xlsx	March 7, 2024, 11:34:35 (UTC+05:45)	27.9 KB	Standard

B. unzipped/ and transformation/ folders were created after trigger

Amazon S3 > Buckets > rrf-bucket

rrf-bucket [Info](#)

Objects Properties Permissions Metrics Management Access Points

Objects (4) [Info](#) Refresh Copy S3 URI Copy URL Download Open Delete Actions Create folder Upload

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

< 1 > ⌕

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	excelfiles/	Folder	-	-	-
<input type="checkbox"/>	transformation/	Folder	-	-	-
<input type="checkbox"/>	unzipped/	Folder	-	-	-
<input type="checkbox"/>	zipfiles/	Folder	-	-	-

Amazon S3> Buckets> rrf-bucket> unzipped/

unzipped/

Copy S3 URI

Objects | Properties

Objects (9) Info

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

Find objects by prefix

< 1 > ⌕

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	RXNATOMARCHIVE.RRF	RRF	March 10, 2024, 16:04:16 (UTC+05:45)	71.4 MB	Standard
<input type="checkbox"/>	RXNCONSO.RRF	RRF	March 10, 2024, 16:04:26 (UTC+05:45)	118.6 MB	Standard
<input type="checkbox"/>	RXNCUI.RRF	RRF	March 10, 2024, 16:04:28 (UTC+05:45)	1.7 MB	Standard
<input type="checkbox"/>	RXNCUICHANGES.RRF	RRF	March 10, 2024, 16:04:27 (UTC+05:45)	14.9 KB	Standard
<input type="checkbox"/>	RXNDOC.RRF	RRF	March 10, 2024, 16:04:28 (UTC+05:45)	214.2 KB	Standard
<input type="checkbox"/>	RXNREL.RRF	RRF	March 10, 2024, 16:05:21 (UTC+05:45)	484.4 MB	Standard
<input type="checkbox"/>	RXNSAB.RRF	RRF	March 10, 2024, 16:05:27 (UTC+05:45)	9.8 KB	Standard
<input type="checkbox"/>	RXNSAT.RRF	RRF	March 10, 2024, 16:06:09 (UTC+05:45)	498.7 MB	Standard
<input type="checkbox"/>	RXNSTY.RRF	RRF	March 10, 2024, 16:06:16 (UTC+05:45)	18.4 MB	Standard

Amazon S3> Buckets> rrf-bucket> transformation/

transformation/

Copy S3 URI

Objects | Properties

Objects (9) Info

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

Find objects by prefix

< 1 > ⌕

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	RXNATOMARCHIVE.txt	txt	March 10, 2024, 16:04:15 (UTC+05:45)	63.2 MB	Standard
<input type="checkbox"/>	RXNCONSO.txt	txt	March 10, 2024, 16:04:25 (UTC+05:45)	118.7 MB	Standard
<input type="checkbox"/>	RXNCUI.txt	txt	March 10, 2024, 16:04:28 (UTC+05:45)	1.6 MB	Standard
<input type="checkbox"/>	RXNCUICHANGES.txt	txt	March 10, 2024, 16:04:27 (UTC+05:45)	15.0 KB	Standard
<input type="checkbox"/>	RXNDOC.txt	txt	March 10, 2024, 16:04:28 (UTC+05:45)	211.0 KB	Standard
<input type="checkbox"/>	RXNREL.txt	txt	March 10, 2024, 16:05:15 (UTC+05:45)	522.3 MB	Standard
<input type="checkbox"/>	RXNSAB.txt	txt	March 10, 2024, 16:05:27 (UTC+05:45)	10.0 KB	Standard
<input type="checkbox"/>	RXNSAT.txt	txt	March 10, 2024, 16:06:03 (UTC+05:45)	497.9 MB	Standard
<input type="checkbox"/>	RXNSTY.txt	txt	March 10, 2024, 16:06:16 (UTC+05:45)	18.1 MB	Standard

C. Here is the RXNATOMARCHIVE.rrf files after transformation of dates in the txt format

Raw

Formatted

RXAU AU STR ARCHIVE_TIMESTAMP CREATED_TIMESTAMP UPDATED_TIMESTAMP CODE IS_BRAND LAT LAST_RELEASED SAUI VSAB RXCU SAB TTY MERGED_TO_RXCU Code Set Version Month
947 A10335796 Mesna 2020-04-27 2005-03-10 2020-04-27 44 ENG 2020-04-06 RXNORM_19AB_200406F 44 RXNORM IN 44 RxNorm 2024-02-05
1424 A10334758 beta-Alanine 2020-04-27 2005-03-10 2020-04-27 61 ENG 2020-04-06 RXNORM_19AB_200406F 61 RXNORM IN 61 RxNorm 2024-02-05
1684 A10334529 4-Aminobenzoic Acid 2020-04-27 2005-03-10 2020-04-27 74 ENG 2020-04-06 RXNORM_19AB_200406F 74 RXNORM IN 74 RxNorm 2024-02-05
2192 A16791816 Eicosapentaenoic Acid 2020-04-27 2005-03-10 2020-04-27 90 ENG 2020-04-06 RXNORM_19AB_200406F 90 RXNORM PIN 90 RxNorm 2024-02-05
2265 A10334531 5-Hydroxytryptophan 2020-04-27 2005-03-10 2020-11-06 94 ENG 2020-04-06 RXNORM_19AB_200406F 94 RXNORM IN 94 RxNorm 2024-02-05
2311 A16793037 Ticlopidine Hydrochloride 2020-04-27 2005-03-10 2020-04-27 97 ENG 2020-04-06 RXNORM_19AB_200406F 97 RXNORM PIN 97 RxNorm 2024-02-05
2332 A10334533 6-Aminocaproic Acid 2020-04-27 2005-03-10 2020-04-27 99 ENG 2020-04-06 RXNORM_19AB_200406F 99 RXNORM IN 99 RxNorm 2024-02-05
2453 A10334534 6-Mercaptopurine 2010-10-21 2005-03-10 2010-10-21 103 ENG 2010-10-04 RXNORM_10AA_101004F 103 RXNORM IN 103 RxNorm 2024-02-05
2663 A10336065 Oxyquinoline 2020-04-27 2005-03-10 2020-04-27 110 ENG 2020-04-06 RXNORM_19AB_200406F 110 RXNORM IN 110 RxNorm 2024-02-05
4330 A10334539 Acebutolol 2020-04-27 2005-03-10 2020-04-27 149 ENG 2020-04-06 RXNORM_19AB_200406F 149 RXNORM IN 149 RxNorm 2024-02-05
4414 A10334540 Acenocoumarol 2020-04-27 2005-03-10 2020-04-27 154 ENG 2020-04-06 RXNORM_19AB_200406F 154 RXNORM IN 154 RxNorm 2024-02-05
4458 A10334542 Acepromazine 2020-04-27 2005-03-10 2020-04-27 155 ENG 2020-04-06 RXNORM_19AB_200406F 155 RXNORM IN 155 RxNorm 2024-02-05
4565 A10334544 Acetanilide 2006-11-15 2005-03-10 2006-11-15 162 ENG 0000-00-00 RXNORM_06AC_061012F 162 RXNORM IN 162 RxNorm 2024-02-05
4655 A10334545 Acetazolamide 2020-04-27 2005-03-10 2020-04-27 167 ENG 2020-04-06 RXNORM_19AB_200406F 167 RXNORM IN 167 RxNorm 2024-02-05
4714 A10338597 Acetic Acid 2020-04-27 2005-03-10 2020-04-27 168 ENG 2020-04-06 RXNORM_19AB_200406F 168 RXNORM IN 168 RxNorm 2024-02-05
4855 A10334546 Acetohexamide 2020-04-27 2005-03-10 2020-04-27 173 ENG 2020-04-06 RXNORM_19AB_200406F 173 RXNORM IN 173 RxNorm 2024-02-05
4998 A10334547 Acetone 2020-04-27 2005-03-10 2020-04-27 178 ENG 2020-04-06 RXNORM_19AB_200406F 178 RXNORM IN 178 RxNorm 2024-02-05

D. Here is the RXNSAB.rrf files transformed after adding SVER column and changing the date format

▶	2024-03-10T17:32:09.994+05:45	Row count before transformation: 7222404
▶	2024-03-10T17:32:09.995+05:45	Version month: 2024-02-05
▶	2024-03-10T17:32:10.096+05:45	Row count of RXNSAT after transformation: 7222404
▶	2024-03-10T17:33:08.069+05:45	The RXNSTY.RRF is read from zip file.
▶	2024-03-10T17:33:08.811+05:45	Row count before transformation: 461874
▶	2024-03-10T17:33:08.812+05:45	Version month: 2024-02-05
▶	2024-03-10T17:33:08.818+05:45	Row count of RXNSTY after transformation: 461874