

## Lambda Structure:

Lambda > Functions > rxnrom\_process

### rxnrom\_process

▼ Function overview Info

Diagram Template

rxnrom\_process

Layers (2)

S3

+ Add trigger

+ Add destination

```
1 import boto3
2 import os
3 import zipfile
4 import tempfile
5 import time
6 import pandas as pd
7 import openpyxl
8 from datetime import datetime
9 import logging
10
11 logger = logging.getLogger()
12 logger.setLevel(logging.INFO)
13
14 s3 = boto3.client('s3')
15
16 def extract_files_from_zip(zip_file_path, destination_path):
17     with zipfile.ZipFile(zip_file_path, 'r') as zip_ref:
18         zip_ref.extractall(destination_path)
19
20 def read_excel_sheets(file_path):
21     try:
22         wb = openpyxl.load_workbook(file_path)
23         return wb.sheetnames
24     except Exception as e:
25         print(f"Error reading Excel file: {e}")
26         return []
27
28 def extract_first_column_values(sheet):
29     column_values = []
30     for row in sheet.iter_rows(min_row=1, max_row=sheet.max_row, min_col=1, max_col=1, values_only=True):
31         column_values.extend(row)
32     return column_values
33
34 #Functions to update the date format
35 def update_RXISAB_df(df):
36     df['VSTART'] = df['VSTART'].apply(lambda x: datetime.strptime(x, '%Y-%m-%d').strftime('%Y-%m-%d') if pd.notnull(x) else x)
37     df['VEND'] = df['VEND'].apply(lambda x: datetime.strptime(x, '%Y-%m-%d').strftime('%Y-%m-%d') if pd.notnull(x) else x)
```

9:15 Python Spaces: 4

```
Environment
  rxnorm_process - /
  lambda_function.py

38 def update_RXNORMARCHIVE_df(df):
39     df['CREATED_TIMESTAMP'] = df['CREATED_TIMESTAMP'].apply(lambda x: datetime.strptime(x, "%m/%d/%Y %I:%M:%S %p").strftime("%Y-%m-%d %I:%M:%S %p") if
40     df['UPDATED_TIMESTAMP'] = df['UPDATED_TIMESTAMP'].apply(lambda x: datetime.strptime(x, "%m/%d/%Y %I:%M:%S %p").strftime("%Y-%m-%d %I:%M:%S %p") if
41     df['LAST_RELEASED'] = df['LAST_RELEASED'].fillna('').astype(str).apply(lambda x: datetime.strptime(x, '%d-%b-%y').strftime('%Y-%m-%d') if x else ''
42
43 # function to set the default value to the columns:
44 def set_default_columns(df):
45     df['Code set'] = 'Rxnorm'
46     df['Version Month'] = version
47
48 # function to extract file name and version of the file
49 def extract_filename_and_version(key):
50     global file_name_parts
51     file_name_parts = key.split("/")[-1].split(".")[0]
52     version_ext = file_name_parts.split("_")[-1]
53     date_object = datetime.strptime(version_ext, "%m%d%Y")
54     global version
55     version = date_object.strftime("%Y-%m-%d")
56     logger.log(logging.INFO, "This is version: %s", version)
57
58
59 # function to unzip the .zip file and save to 'unzipped' folder
60 def unzip_file_saveto_unzipped(key):
61     # Check if the uploaded file is a zip file
62     if key.endswith('.zip'):
63         # Create a temporary directory to extract the zip file contents
64         with tempfile.TemporaryDirectory() as tmp_dir:
65             # Download the zip file to the temporary directory
66             local_zip_file = os.path.join(tmp_dir, 'upload.zip')
67             s3.download_file(bucket, key, local_zip_file)
68
69             # Extract the contents of the zip file
70             extract_files_from_zip(local_zip_file, tmp_dir)
```

```
Environment
  rxnorm_process - /
  lambda_function.py

75     extracted_file = s3.get_object(Bucket=bucket, Key=key)
76     extracted_file_path = os.path.join(rrf_folder, extracted_file)
77     unzipped_key = f'unzipped/{extracted_file}'
78     s3.upload_file(extracted_file_path, bucket, unzipped_key)
79     print(f"Successfully extracted and stored contents of {key} in 'unzipped' folder.")
80 else:
81     print(f"No 'rrf' folder found in the extracted files from {key}.")
82
83
84 # function to work with excel file:
85 def excel_read(bucket):
86     keye = 'upload/RxNorm_Header.xlsx'
87     # Download the Excel file to the /tmp directory
88     local_excel_file = '/tmp/RxNorm_Header.xlsx'
89     s3.download_file(bucket, keye, local_excel_file)
90
91     # Check if the Excel file exists
92     if os.path.exists(local_excel_file):
93
94         try:
95             # Read all the sheets of the Excel file
96             sheet_names = read_excel_sheets(local_excel_file)
97             print(f"Sheets in {keye}:")
98             for sheet_name in sheet_names:
99                 # Print(sheet_name)
100
101             # Load the sheet
102             wb = openpyxl.load_workbook(local_excel_file)
103             sheet = wb[sheet_name]
104
105             # Extract first column values
106             first_column_values = extract_first_column_values(sheet)
107
108             # Create a list with name = 'sheet_name-Headers'
109             list_name = f'{sheet_name}.Headers'
110             globals()[list_name] = first_column_values
111
112             print(f"List {list_name} created with values: {first_column_values}")
```

9:15 Python Spaces: 4

Environment

rxnorm\_process

lambda\_function.py

lambda\_function.py

```

120 # function to read the files from the unzipped folder, transform and save to destination
121 def read_transform_save():
122     unzipped_folder_prefix = 'unzipped/'
123     global destination_folder
124     destination_folder = 'destination/'
125
126     # List objects in the 'unzipped' folder
127     response = s3.list_objects_v2(Bucket=bucket, Prefix=unzipped_folder_prefix)
128
129     # Check if the 'unzipped' folder exists
130     if 'Contents' in response:
131         for obj in response['Contents']:
132             object_key = obj['Key']
133
134             # Extract the file name from the object key
135             file_name = os.path.basename(object_key)
136
137             # Check if the file is an .RRF file
138             if file_name.endswith('.RRF'):
139                 try:
140                     # Get the object from S3
141                     response = s3.get_object(Bucket=bucket, Key=object_key)
142
143                     # Read the file into a DataFrame without header
144                     df = pd.read_csv(response['Body'], sep='|', header=None)
145                     # deleting the extra null column
146                     df = df.iloc[:, :-1]
147                     set_default_columns(df)
148                     # count the records before transformation
149                     global record_count_before
150                     record_count_before = df.shape[0]
151                     # Add headers to the pandas DataFrame
152                     list_name = f"{os.path.splitext(file_name)[0]}_Headers"
153                     df.columns = globals()[list_name]

```

Go to Anything (Ctrl-P)

Environment

rxnorm\_process

lambda\_function.py

lambda\_function.py

lambda\_function.py

Environment Var

```

162 update_RXNNAI(OPAHLHIVE_df(df)
163 #calling function to convert the df to csv and save to destination folder
164 convert_save(df, file_name)
165
166 except Exception as e:
167     print(f"Error reading or uploading file {file_name}: {e}")
168 else:
169     print(f"No files found in the 'unzipped' folder.")
170
171
172 # function to convert df to csv and save to destination folder
173 def convert_save(df, file_name):
174     record_count_after = df.shape[0]
175     # Convert DataFrame to CSV with comma delimiter
176     csv_buffer = df.to_csv(index=False, sep=',')
177
178     # Upload CSV file to S3 with the same name as the original file
179     s3.put_object(Bucket=bucket, Key=destination_folder + file_name.replace('.RRF', '.csv'), Body=csv_buffer)
180
181     print(f"CSV file '{file_name.replace('.RRF', '.csv')}' uploaded to '{destination_folder}' in S3.
182           Records of {file_name} before transformation: {record_count_before}
183           Records of {file_name} after transformation: {record_count_after}")
184
185 # main trigger handler
186 def lambda_handler(event, context):
187     for record in event['Records']:
188         global bucket
189         bucket = record['s3']['bucket']['name']
190         key = record['s3']['object']['key']
191         extract_filename_and_version(key)
192         # unzip the file and store to 'unzipped' folder
193         unzip_file_saveto_unzipped(key)
194         #working with excel file and read the sheets
195         excel_read(bucket)
196         # Read files from the 'unzipped' folder transform it and save to destination
197         read_transform_save()

```

## S3 folder structure:

The screenshot shows the Amazon S3 console interface for the bucket **rxnorm1313**. The breadcrumb navigation at the top indicates the path: [Amazon S3](#) > [Buckets](#) > [rxnorm1313](#). Below the bucket name, there are tabs for **Objects**, **Properties**, **Permissions**, **Metrics**, **Management**, and **Access Points**. The **Objects** tab is active, displaying a list of objects. Above the list, there are buttons for **Refresh**, **Copy S3 URI**, **Copy URL**, **Download**, **Open**, **Delete**, and **Actions**. A search bar with the placeholder text "Find objects by prefix" is present. The table below lists three folders: **destination/**, **unzipped/**, and **upload/**, all of which are of type **Folder** and have a **Last modified** date of **-** and a **Size** of **-**.

| <input type="checkbox"/> | Name                         | Type   | Last modified | Size |
|--------------------------|------------------------------|--------|---------------|------|
| <input type="checkbox"/> | <a href="#">destination/</a> | Folder | -             | -    |
| <input type="checkbox"/> | <a href="#">unzipped/</a>    | Folder | -             | -    |
| <input type="checkbox"/> | <a href="#">upload/</a>      | Folder | -             | -    |

## Upload Folder:

The screenshot shows the Amazon S3 console interface for the **upload/** folder within the **rxnorm1313** bucket. The breadcrumb navigation at the top indicates the path: [Amazon S3](#) > [Buckets](#) > [rxnorm1313](#) > [upload/](#). Below the folder name, there are tabs for **Objects** and **Properties**. The **Objects** tab is active, displaying a list of objects. Above the list, there are buttons for **Refresh** and **Actions**. A search bar with the placeholder text "Find objects by prefix" is present. The table below lists two objects: **RxNorm\_full\_02052024.zip** and **RxNorm\_Header.xlsx**, both of which are of type **zip** and **xlsx** respectively and have a **Last modified** date of **M**.

| <input type="checkbox"/> | Name                                     | Type | Last modified |
|--------------------------|------------------------------------------|------|---------------|
| <input type="checkbox"/> | <a href="#">RxNorm_full_02052024.zip</a> | zip  | M             |
| <input type="checkbox"/> | <a href="#">RxNorm_Header.xlsx</a>       | xlsx | M             |

Unzipped folder:

Amazon S3

>

Buckets

>

rxnorm1313

>

unzipped/

unzipped/

Objects

Properties

Objects (9)

Info

Copy S3 URI

Copy URL

Download

Open

Delete

Actions

Create folder

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 Inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

Find objects by prefix

| <input type="checkbox"/> | Name                               | Type | Last modified                       | Size     | Storage class |
|--------------------------|------------------------------------|------|-------------------------------------|----------|---------------|
| <input type="checkbox"/> | <a href="#">RXNATOMARCHIVE.RRF</a> | RRF  | March 7, 2024, 00:48:55 (UTC+05:45) | 71.4 MB  | Standard      |
| <input type="checkbox"/> | <a href="#">RXNCONSO.RRF</a>       | RRF  | March 7, 2024, 00:48:56 (UTC+05:45) | 118.6 MB | Standard      |
| <input type="checkbox"/> | <a href="#">RXNCUI.RRF</a>         | RRF  | March 7, 2024, 00:48:55 (UTC+05:45) | 1.7 MB   | Standard      |
| <input type="checkbox"/> | <a href="#">RXNCUICHANGES.RRF</a>  | RRF  | March 7, 2024, 00:48:55 (UTC+05:45) | 14.9 KB  | Standard      |
| <input type="checkbox"/> | <a href="#">RXNDOC.RRF</a>         | RRF  | March 7, 2024, 00:48:55 (UTC+05:45) | 214.2 KB | Standard      |
| <input type="checkbox"/> | <a href="#">RXNREL.RRF</a>         | RRF  | March 7, 2024, 00:48:58 (UTC+05:45) | 484.4 MB | Standard      |
| <input type="checkbox"/> | <a href="#">RXNSAB.RRF</a>         | RRF  | March 7, 2024, 00:48:55 (UTC+05:45) | 9.8 KB   | Standard      |
| <input type="checkbox"/> | <a href="#">RXNSAT.RRF</a>         | RRF  | March 7, 2024, 00:48:50 (UTC+05:45) | 498.7 MB | Standard      |
| <input type="checkbox"/> | <a href="#">RXNSTY.RRF</a>         | RRF  | March 7, 2024, 00:48:50 (UTC+05:45) | 18.4 MB  | Standard      |

Destination folder:

Amazon S3

>

Buckets

>

rxnorm1313

>

destination/

destination/

Objects

Properties

Objects (9)

Info

Copy S3 URI

Copy URL

Download

Open

Delete

Actions

Create folder

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 Inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

Find objects by prefix

| <input type="checkbox"/> | Name                               | Type | Last modified                       | Size     | Storage class |
|--------------------------|------------------------------------|------|-------------------------------------|----------|---------------|
| <input type="checkbox"/> | <a href="#">RXNATOMARCHIVE.csv</a> | csv  | March 7, 2024, 00:49:21 (UTC+05:45) | 77.8 MB  | Standard      |
| <input type="checkbox"/> | <a href="#">RXNCONSO.csv</a>       | csv  | March 7, 2024, 00:49:34 (UTC+05:45) | 138.5 MB | Standard      |
| <input type="checkbox"/> | <a href="#">RXNCUI.csv</a>         | csv  | March 7, 2024, 00:49:35 (UTC+05:45) | 2.1 MB   | Standard      |
| <input type="checkbox"/> | <a href="#">RXNCUICHANGES.csv</a>  | csv  | March 7, 2024, 00:49:35 (UTC+05:45) | 17.9 KB  | Standard      |
| <input type="checkbox"/> | <a href="#">RXNDOC.csv</a>         | csv  | March 7, 2024, 00:49:35 (UTC+05:45) | 271.7 KB | Standard      |
| <input type="checkbox"/> | <a href="#">RXNREL.csv</a>         | csv  | March 7, 2024, 00:50:43 (UTC+05:45) | 645.2 MB | Standard      |
| <input type="checkbox"/> | <a href="#">RXNSAB.csv</a>         | csv  | March 7, 2024, 00:50:50 (UTC+05:45) | 10.3 KB  | Standard      |
| <input type="checkbox"/> | <a href="#">RXNSAT.csv</a>         | csv  | March 7, 2024, 00:51:53 (UTC+05:45) | 622.4 MB | Standard      |
| <input type="checkbox"/> | <a href="#">RXNSTY.csv</a>         | csv  | March 7, 2024, 00:52:02 (UTC+05:45) | 26.1 MB  | Standard      |

Visualizing the original file and Transformed file  
Example for RXNSAB file,

RXNSAB.RRF:

Query results

Query results are not available after you choose Close or navigate away. Choose Download results to download a copy of the following query results.

Download results

Status

Successfully returned 13 records in 2149 ms

Bytes returned: 10066 B

Raw

Formatted

|          |          |                        |             |                                                       |           |              |        |                                                          |                                                                       |                              |          |        |
|----------|----------|------------------------|-------------|-------------------------------------------------------|-----------|--------------|--------|----------------------------------------------------------|-----------------------------------------------------------------------|------------------------------|----------|--------|
| C5233827 | C1140218 | MMSL_2024_01_01        | MMSL        | Multum MediSource Lexicon                             | MMSL      | 2024_01_01   | 2020AA | Multum Information Services                              | 3200 Cherry Creek South Drive                                         |                              |          |        |
| C5233830 | C1140182 | MMX_2024_01_02         | MMX         | Micromedex RED BOOK                                   | MMX       | 2024_01_02   | 2020AA | Micromedex                                               | 6200 South Syracuse Way                                               |                              |          |        |
| C5233828 | C1140261 | NDDF_2024_01_03        | NDDF        | FDB MedKnowledge (formerly NDDF Plus)                 | NDDF      | 2024_01_03   | 2020AA | First Databank Customer Support                          | 701 Gateway Blvd                                                      |                              |          |        |
| C5233835 | C1140284 | RXNORM_20AA_240205F    | RXNORM      | RxNorm Vocabulary                                     | RXNORM    | 20AA_240205F | 2020AA | RxNorm Customer Service                                  | U.S. National Library of Medicine                                     | 8600 Rockville Pike          | Bethesda | MD     |
| C5233831 | C2720507 | SNOMEDCT_US_2023_06_30 | SNOMEDCT_US | US Edition of SNOMED CT                               | SNOMEDCT  | 2023_06_30   | 2020AA | National Library Of Medicine                             | NLM is a Charter Member of SNOMED International on behalf of the U.S. | National Library of Medicine | 86       |        |
| C5233837 | C1140288 | VANDF_2023_12_29       | VANDF       | Veterans Health Administration National Drug File     | VANDF     | 2023_12_29   | 2020AA | Michael Lincoln                                          |                                                                       |                              |          |        |
| C5233836 | C1812643 | MTHSPL_2024_01_27      | MTHSPL      | Metathesaurus FDA Structured Product Labels           | MTHSPL    | 2024_01_27   | 2020AA | RxNorm Customer Service                                  | U.S. National Library of Medicine                                     | 8600 Rockville Pike          | Bethesda | MD     |
| C5233833 | C1876239 | GS_2024_01_05          | GS          | Gold Standard Drug Database                           | GS        | 2024_01_05   | 2020AA | RxNorm Customer Service                                  | U.S. National Library of Medicine                                     | 8600 Rockville Pike          | Bethesda | MD     |
| C5233831 | C4722517 | ATC_2024_01_19         | ATC         | Anatomical Therapeutic Chemical Classification System | ATC       | 2024_01_19   | 2020AA | WHO Collaborating Centre for Drug Statistics Methodology | Norwegian Institute of Public Health                                  | P.O.Box 4404 Nydalen         | Oslo     | Norway |
| C5233832 | C3539983 | CVX_2024_01_08         | CVX         | Vaccines Administered                                 | CVX       | 2024_01_08   | 2020AA | CDC                                                      |                                                                       |                              |          |        |
| C5233829 | C3858951 | MTHCMSFRF_2020         | MTHCMSFRF   | Metathesaurus CMS Formulary Reference File            | MTHCMSFRF | 2020         | 2020AA | RxNorm Customer Service                                  | U.S. National Library of Medicine                                     | 8600 Rockville Pike          | Bethesda | MD     |
| C5233834 | C4554231 | USP_2024_01_08         | USP         | USP Compndial Nomenclature                            | USP       | 2024_01_08   | 2020AA | Jeffrey Shi                                              |                                                                       |                              |          |        |

RXNSAB.csv:

bytes returned: 10066 B

Raw

Formatted

|          |          |                        |             |                                                       |           |                |            |      |        |       |             |
|----------|----------|------------------------|-------------|-------------------------------------------------------|-----------|----------------|------------|------|--------|-------|-------------|
| VCUI     | RCUI     | VSAB                   | RSAB        | SON                                                   | SF        | SVER           | VSTART     | VEND | IMETA  | RMETA | SLC         |
| C5233827 | C1140218 | MMSL_2024_01_01        | MMSL        | Multum MediSource Lexicon                             | MMSL      | 2024_01_01     |            |      | 2020AA |       | Multum In   |
| C5233830 | C1140182 | MMX_2024_01_02         | MMX         | Micromedex RED BOOK                                   | MMX       | 2024_01_02     |            |      | 2020AA |       | Micromed    |
| C5233828 | C1140261 | NDDF_2024_01_03        | NDDF        | FDB MedKnowledge (formerly NDDF Plus)                 | NDDF      | 2024_01_03     |            |      | 2020AA |       | First Data  |
| C5233835 | C1140284 | RXNORM_20AA_240205F    | RXNORM      | RxNorm Vocabulary                                     | RXNORM    | 20AA_240205F   |            |      | 2020AA |       | RxNorm Cl   |
| C5233831 | C2720507 | SNOMEDCT_US_2023_06_30 | SNOMEDCT_US | US Edition of SNOMED CT                               | SNOMEDCT  | 2023_06_30     |            |      | 2020AA |       | National Li |
| C5233837 | C1140288 | VANDF_2023_12_29       | VANDF       | Veterans Health Administration National Drug File     | VANDF     | 2023_12_29     |            |      | 2020AA |       | Michael Lir |
| C5233836 | C1812643 | MTHSPL_2024_01_27      | MTHSPL      | Metathesaurus FDA Structured Product Labels           | MTHSPL    | 2024_01_27     | 2006-12-21 |      | 2020AA |       | RxNorm Cl   |
| C5233833 | C1876239 | GS_2024_01_05          | GS          | Gold Standard Drug Database                           | GS        | 2024_01_05     | 2007-05-04 |      | 2020AA |       | RxNorm Cl   |
| C5233831 | C4722517 | ATC_2024_01_19         | ATC         | Anatomical Therapeutic Chemical Classification System | ATC       | 2024_01_19     |            |      | 2020AA |       | WHO Coll    |
| C5233832 | C3539983 | CVX_2024_01_08         | CVX         | Vaccines Administered                                 | CVX       | 2024_01_08     |            |      | 2020AA |       | CDC, Nati   |
| C5233829 | C3858951 | MTHCMSFRF_2020         | MTHCMSFRF   | Metathesaurus CMS Formulary Reference File            | MTHCMSFRF | 2020           |            |      | 2020AA |       | RxNorm Cl   |
| C5233838 | C4255544 | DRUGBANKS_0_2024_01_04 | DRUGBANK    | DrugBank                                              | DRUGBANK  | 5_0_2024_01_04 |            |      | 2020AA |       | OMx Pers    |
| C5233834 | C4554231 | USP_2024_01_08         | USP         | USP Compndial Nomenclature                            | USP       | 2024_01_08     |            |      | 2020AA |       | Jeffrey Shi |

Validating the record for RXNSAB file:

|   |                               |                                                         |
|---|-------------------------------|---------------------------------------------------------|
| ▶ | 2024-03-07T00:50:49.064+05:45 | CSV file 'RXNSAB.csv' uploaded to 'destination/' in S3. |
| ▶ | 2024-03-07T00:50:49.064+05:45 | Records of RXNSAB.RRF before transformation: 13         |
| ▶ | 2024-03-07T00:50:49.064+05:45 | Records of RXNSAB.RRF after transformation: 13          |

Same goes for all the 9 files which can be verified with cloudwatch logs.