

# RxNorm Lambda project

## General overview

### 1. Create the lambda function with s3 bucket trigger

[Lambda](#) > [Functions](#) > RxNorm\_Lambda\_test

### RxNorm\_Lambda\_test

Throttle Copy ARN Actions

Function overview Info

Export to Application Composer Download

Diagram Template

RxNorm\_Lambda\_test

Layers (1)

S3

+ Add trigger

+ Add destination

Description

-

Last modified

1 hour ago

Function ARN

arn:aws:lambda:us-east-1:866388144037:function:RxNorm\_Lambda\_test

Function URL Info

### 2. Add necessary configuration like memory, timeout and layers for lambda

Configuration Aliases Versions

### General configuration Info

Edit

Description	Memory	Ephemeral storage
-	6024 MB	6024 MB
Timeout	SnapStart Info	
8 min 3 sec	None	

### Layers Info

Edit Add a layer

Merge order	Name	Layer version	Compatible runtimes	Compatible architectures	Version ARN
1	AWSSDKPandas-Python312	4	python3.12	x86_64	arn:aws:lambda:us-east-1:336392948345:layer:AWSSDKPandas-Python312:4

### 3. This is my final working structure

[Amazon S3](#) > [Buckets](#) > rxnorm-s3bucket-test

### rxnorm-s3bucket-test Info

Objects Properties Permissions Metrics Management Access Points

Objects (4) Info


Copy Copy S3 URI Copy URL Download Open Delete Actions Create folder Upload

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 Inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)


Find objects by prefix

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	excel_file/	Folder	-	-	-
<input type="checkbox"/>	extracted_files/	Folder	-	-	-
<input type="checkbox"/>	transformed_files/	Folder	-	-	-
<input type="checkbox"/>	zipfile/	Folder	-	-	-

4. First upload the [RxNorm\\_full\\_02052024.zip](#) in the zip file which will trigger the lambda and executes the specific code



5. Upload the [RxNorm\\_Header.xlsx](#) in excel\_file folder to read the headers



```
def read_excel_from_s3(bucket):
    try:
        folder_path = 'excel_file/'
        excel_file_name = 'RxNorm_Header.xlsx'
        key = folder_path + excel_file_name

        # Download the Excel file to the /tmp directory
        local_excel_file = '/tmp/RxNorm_Header.xlsx'

        s3.download_file(bucket, key, local_excel_file)

        # Check if the Excel file exists
        if os.path.exists(local_excel_file):
            print(f"Excel file downloaded to: {local_excel_file}")
        else:
            print(f"Excel file is not found yet")

        # Read the Excel file into an ExcelFile object
        excel_file = pd.ExcelFile(local_excel_file)
        print(f"The excel file RxNorm_Header.xlsx is read from S3")

        # Get the sheet names
        sheet_names = excel_file.sheet_names
        headers_data = {}
        for sheets in sheet_names:
            print(sheets)
            sheets_data = excel_file.parse(sheets, header=None)
            excel_headers = sheets_data.iloc[:, 0].tolist()
            print(f"The Sheet, {sheets} is read from RxNorm_Header.xlsx, for it's header {excel_headers}")
            headers_data[sheets] = excel_headers

        return headers_data
    except Exception as e:
```

Code that reads the headers for each file when the lambda triggers.

6. Now when i upload the zip file it will only extract the file in the RRF folder and dumps in the extracted\_files folder

Amazon S3 > Buckets > rxnorm-s3bucket-test > extracted\_files/

extracted\_files/ Copy S3 URI

Objects | Properties

Objects (9) [Info](#) Copy S3 URI Copy URL Download Open Delete Actions Create folder Upload

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 Inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	<a href="#">RXNATOMARCHIVE.RRF</a>	RRF	March 11, 2024, 15:03:56 (UTC+05:45)	71.4 MB	Standard
<input type="checkbox"/>	<a href="#">RXNCONSO.RRF</a>	RRF	March 11, 2024, 15:03:57 (UTC+05:45)	118.6 MB	Standard
<input type="checkbox"/>	<a href="#">RXNCUI.RRF</a>	RRF	March 11, 2024, 15:03:58 (UTC+05:45)	1.7 MB	Standard
<input type="checkbox"/>	<a href="#">RXNCUICHANGES.RRF</a>	RRF	March 11, 2024, 15:03:58 (UTC+05:45)	14.9 KB	Standard
<input type="checkbox"/>	<a href="#">RXNDOC.RRF</a>	RRF	March 11, 2024, 15:03:58 (UTC+05:45)	214.2 KB	Standard
<input type="checkbox"/>	<a href="#">RXNREL.RRF</a>	RRF	March 11, 2024, 15:04:00 (UTC+05:45)	484.4 MB	Standard
<input type="checkbox"/>	<a href="#">RXNSAB.RRF</a>	RRF	March 11, 2024, 15:04:04 (UTC+05:45)	9.8 KB	Standard
<input type="checkbox"/>	<a href="#">RXNSAT.RRF</a>	RRF	March 11, 2024, 15:04:07 (UTC+05:45)	498.7 MB	Standard
<input type="checkbox"/>	<a href="#">RXNSTY.RRF</a>	RRF	March 11, 2024, 15:04:12 (UTC+05:45)	18.4 MB	Standard

```
#Read the zip file and relocate the files in the rrf folder
def read_and_relocate_rrf_files(s3, bucket, key):
    try:
        # My columns for the Version Month field
        text_data = key
        text_file = text_data.split('/')[1].split('.')[2].split('.')[0]
        date_string = datetime.strptime(text_file, "%m%d%Y")
        date_object = date_string.strftime('%B %d, %Y')
        version_month.append(date_object)
        print(f"The data for Version Month is: {version_month[0]}")

        # Reading the rrf files
        zip_response = s3.get_object(Bucket=bucket, Key=key)
        zip_data = zip_response['Body'].read()
        file_path = 'rrf'
        destination_path = 'extracted_files'

        # Read the zip file without unzipping it
        with zipfile.ZipFile(io.BytesIO(zip_data)) as zip_ref:
            for file_info in zip_ref.infolist():
                if file_info.filename.startswith(file_path) and not file_info.filename.endswith('/'):
                    filename = os.path.basename(file_info.filename)
                    if filename.endswith('.RRF'):
                        print(f"The {filename} is read from zip file.")
                        with zip_ref.open(file_info) as source_file:
                            s3.put_object(
                                Bucket=bucket,
                                Key=f"{destination_path}/{filename}",
                                Body=source_file.read()
                            )
                        print(f"The {filename} is relocated to {destination_path}")

        return {
            'statusCode': 200,
            'body': json.dumps('Zip file processed successfully')
        }
    except Exception as e:
```

This code will also calculate the data for Version Month column along with relocating the files in the rrf folder to extracted\_file

## 7. Now the below code will update the columns for the required files

```
# Updating the dataframe for some files
def update_dataframe(df, file_name):
    if file_name == 'RXNSAB.RRF':
        print(f'Updating the columns VSTART and VEND for {file_name}....')
        df['VSTART'] = df['VSTART'].apply(lambda x: datetime.strptime(x, '%Y_%m_%d').strftime('%Y-%m-%d') if pd.notnull(x) else x)
        df['VEND'] = df['VEND'].apply(lambda x: datetime.strptime(x, '%Y_%m_%d').strftime('%Y-%m-%d') if pd.notnull(x) else x)
    elif file_name == 'RXNATOMARCHIVE.RRF':
        print(f'Updating the columns CREATED_TIMESTAMP, UPDATED_TIMESTAMP and LAST_RELEASED for {file_name}....')
        df['CREATED_TIMESTAMP'] = df['CREATED_TIMESTAMP'].apply(lambda x: datetime.strptime(x, '%m/%d/%Y %I:%M:%S %p').strftime('%Y-%m-%d %I:%M:%S %p') if pd
        df['UPDATED_TIMESTAMP'] = df['UPDATED_TIMESTAMP'].apply(lambda x: datetime.strptime(x, '%m/%d/%Y %I:%M:%S %p').strftime('%Y-%m-%d %I:%M:%S %p') if pd
        df['LAST_RELEASED'] = df['LAST_RELEASED'].fillna('').astype(str).apply(lambda x: datetime.strptime(x, '%d-%b-%y').strftime('%Y-%m-%d') if x else '')
    return df

# Function for reading the unzipped rrf files and converting it to csv after adding headers
```

## 8. Now read the files from the extracted\_file folder and make necessary transformation and dump into the transform\_file folder

Amazon S3 > Buckets > rxnorm-s3bucket-test > transformed\_files/

transformed\_files/ Copy S3 URI

Objects Properties

Objects (8) Info Refresh Copy S3 URI Copy URL Download Open Delete Actions Create folder Upload

Objects are the fundamental entities stored in Amazon S3. You can use [Amazon S3 Inventory](#) to get a list of all objects in your bucket. For others to access your objects, you'll need to explicitly grant them permissions. [Learn more](#)

Find objects by prefix

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	<a href="#">RXNATOMARCHIVE.CSV</a>	CSV	March 11, 2024, 15:04:26 (UTC+05:45)	81.0 MB	Standard
<input type="checkbox"/>	<a href="#">RXNCONSO.CSV</a>	CSV	March 11, 2024, 15:04:37 (UTC+05:45)	148.2 MB	Standard
<input type="checkbox"/>	<a href="#">RXNCUI.CSV</a>	CSV	March 11, 2024, 15:04:39 (UTC+05:45)	2.4 MB	Standard
<input type="checkbox"/>	<a href="#">RXNCUICHANGES.CSV</a>	CSV	March 11, 2024, 15:04:39 (UTC+05:45)	19.2 KB	Standard
<input type="checkbox"/>	<a href="#">RXNDOC.CSV</a>	CSV	March 11, 2024, 15:04:39 (UTC+05:45)	302.0 KB	Standard
<input type="checkbox"/>	<a href="#">RXNREL.CSV</a>	CSV	March 11, 2024, 15:05:30 (UTC+05:45)	706.6 MB	Standard
<input type="checkbox"/>	<a href="#">RXNSAB.CSV</a>	CSV	March 11, 2024, 15:05:43 (UTC+05:45)	10.4 KB	Standard
<input type="checkbox"/>	<a href="#">RXNSAT.CSV</a>	CSV	March 11, 2024, 15:06:26 (UTC+05:45)	684.4 MB	Standard

```
# Function for reading the unzipped rrf files and converting it to csv after adding headers
def process_files(bucket, headers_data):
    unzipped_folder = 'extracted_files/'
    destination_folder = 'transformed_files'
    headers = read_excel_from_s3(bucket)

    unzipped_response = s3.list_objects_v2(Bucket=bucket, Prefix=unzipped_folder)

    if 'Contents' in unzipped_response:
        for obj in unzipped_response['Contents']:
            object_key = obj['Key']
            file_name = os.path.basename(object_key)
            print(f"Processing file: {file_name}")
            new_file = file_name.split('.')[0]

            if file_name.endswith('.RRF'):
                try:
                    base_filename = file_name.split('.')[0]
                    file_response = s3.get_object(Bucket=bucket, Key=object_key)
                    df = pd.read_csv(file_response['Body'], sep='|', header=None, low_memory=False)
                    df = df.iloc[:, :-1]
                    df['Code set'] = 'Rxnorm'
                    df['Version Month'] = version_month[0]

                    if base_filename in headers:
                        df.columns = headers[base_filename]
                        df = update_dataframe(df, file_name)

                    before_transformation = df.shape[0]

                    print(f"The count of {file_name} before transformation is: {before_transformation}")

                    csv_buffer = df.to_csv(sep=',', index=False)
```

```

149
150
151         s3.put_object(
152             Bucket=bucket,
153             Key=f"{destination_folder}/{new_file}.CSV",
154             Body=csv_buffer
155         )
156
157
158         print(f"The {file_name} is transformed into CSV and stored into {destination_folder} in S3")
159
160         csv_data = StringIO(csv_buffer)
161         csv_data_length = len(csv_data.readlines()) - 1
162
163         print(f"The count of {file_name} after transformation is: {csv_data_length}")
164
165     except Exception as e:
166         print(f"Error reading file {file_name}: {e}")
167
168     else:
169         print("Zipped file is not found")
170
171
172 def lambda_handler(event, context):
173     bucket = event['Records'][0]['s3']['bucket']['name']
174     key = event['Records'][0]['s3']['object']['key']
175
176     read_and_relocate_rrf_files(s3,bucket,key)
177
178     headers_data = read_excel_from_s3(bucket)
179
180     if headers_data:
181         process_files(bucket,headers_data)
182
183

```

## 9. Now open the transformed file and check the data

Amazon S3 > Buckets > rxnorm-s3bucket-test > transformed\_files/ > RXNATOMARCHIVE.CSV

RXNATOMARCHIVE.CSV Info

[Copy S3 URI](#) [Download](#) [Open](#) [Object actions](#)

**Properties** Permissions Versions

**Object overview**

Owner  
aws-labs-cow697479011703159794

AWS Region  
US East (N. Virginia) us-east-1

Last modified  
March 11, 2024, 15:04:26 (UTC+05:45)

Size  
81.0 MB

Type  
CSV

S3 URI  
s3://rxnorm-s3bucket-test/transformed\_files/RXNATOMARCHIVE.CSV

Amazon Resource Name (ARN)  
arn:aws:s3:::rxnorm-s3bucket-test/transformed\_files/RXNATOMARCHIVE.CSV

Entity tag (ETag)  
bdf8a778b6d34a2569102b38a6abfe0f

Object URL  
[https://rxnorm-s3bucket-test.s3.amazonaws.com/transformed\\_files/RXNATOMARCHIVE.CSV](https://rxnorm-s3bucket-test.s3.amazonaws.com/transformed_files/RXNATOMARCHIVE.CSV)

**Object actions**

- Download as
- Share with a presigned URL
- Calculate total size
- Copy
- Move
- Initiate restore
- Query with S3 Select
- Edit actions
- Rename object
- Edit storage class
- Edit server-side encryption
- Edit metadata
- Edit tags
- Make public using ACL

### Query results

Query results are not available after you choose **Close** or navigate away. Choose **Download results** to download a copy of the following query results.

[Download results](#)

#### Status

Successfully returned 200 records in 2676 ms

Bytes returned: 34638 B

Raw

**Formatted**

RXAU1	AU1	STR	ARCHIVE_TIMESTAMP	CREATED_TIMESTAMP	UPDATED_TIMESTAMP	CODE	IS_BRAND	LAT	LAST_RELEASED	SAU1	VSAB
947	A10335796	Mesna	2020-04-27	2005-03-10 02:03:47 PM	2020-04-27 09:04:10 PM	44		ENG	2020-04-06		RXNORM_19AB_2004
1424	A10334758	beta-Alanine	2020-04-27	2005-03-10 02:03:47 PM	2020-04-27 09:04:16 PM	61		ENG	2020-04-06		RXNORM_19AB_2004
1684	A10334529	4-Aminobenzoic Acid	2020-04-27	2005-03-10 02:03:47 PM	2020-04-27 09:04:14 PM	74		ENG	2020-04-06		RXNORM_19AB_2004
2192	A16791816	Elcosapentaenolic Acid	2020-04-27	2005-03-10 02:03:47 PM	2020-04-27 09:04:18 PM	90		ENG	2020-04-06		RXNORM_19AB_2004
2265	A10334531	5-Hydroxytryptophan	2020-04-27	2005-03-10 02:03:47 PM	2020-11-06 06:11:15 AM	94		ENG	2020-04-06		RXNORM_19AB_2004
2311	A16793037	Ticlopidine Hydrochloride	2020-04-27	2005-03-10 02:03:47 PM	2020-04-27 09:04:19 PM	97		ENG	2020-04-06		RXNORM_19AB_2004
2332	A10334533	6-Aminocaproic Acid	2020-04-27	2005-03-10 02:03:47 PM	2020-04-27 09:04:09 PM	99		ENG	2020-04-06		RXNORM_19AB_2004
2453	A10334534	6-Mercaptopurine	2010-10-21	2005-03-10 02:03:47 PM	2010-10-21 02:10:12 AM	103		ENG	2010-10-04		RXNORM_10AA_1010
2663	A10336065	Oxyquinoline	2020-04-27	2005-03-10 02:03:47 PM	2020-04-27 09:04:09 PM	110		ENG	2020-04-06		RXNORM_19AB_2004
4330	A10334539	Acebutolol	2020-04-27	2005-03-10 02:03:47 PM	2020-04-27 09:04:14 PM	149		ENG	2020-04-06		RXNORM_19AB_2004

														< 1 2 3 4 >	
MESTAMP	UPDATED_TIMESTAMP	CODE	IS_BRAND	LAT	LAST_RELEASED	SAUI	VSAB	RXCUI	SAB	TTY	MERGED_TO_RXCUI	Code set	Version Month		
02:03:47 PM	2020-04-27 09:04:10 PM	44		ENG	2020-04-06		RXNORM_19AB_200406F	44	RXNORM	IN	44	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:16 PM	61		ENG	2020-04-06		RXNORM_19AB_200406F	61	RXNORM	IN	61	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:14 PM	74		ENG	2020-04-06		RXNORM_19AB_200406F	74	RXNORM	IN	74	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:18 PM	90		ENG	2020-04-06		RXNORM_19AB_200406F	90	RXNORM	PIN	90	Rxnorm	February 05, 2024		
02:03:47 PM	2020-11-06 06:11:15 AM	94		ENG	2020-04-06		RXNORM_19AB_200406F	94	RXNORM	IN	94	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:19 PM	97		ENG	2020-04-06		RXNORM_19AB_200406F	97	RXNORM	PIN	97	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:09 PM	99		ENG	2020-04-06		RXNORM_19AB_200406F	99	RXNORM	IN	99	Rxnorm	February 05, 2024		
02:03:47 PM	2010-10-21 02:10:12 AM	103		ENG	2010-10-04		RXNORM_10AA_101004F	103	RXNORM	IN	103	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:09 PM	110		ENG	2020-04-06		RXNORM_19AB_200406F	110	RXNORM	IN	110	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:14 PM	149		ENG	2020-04-06		RXNORM_19AB_200406F	149	RXNORM	IN	149	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:16 PM	154		ENG	2020-04-06		RXNORM_19AB_200406F	154	RXNORM	IN	154	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:11 PM	155		ENG	2020-04-06		RXNORM_19AB_200406F	155	RXNORM	IN	155	Rxnorm	February 05, 2024		
02:03:47 PM	2006-11-15 11:11:45 AM	162		ENG			RXNORM_06AC_061012F	162	RXNORM	IN	162	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:15 PM	167		ENG	2020-04-06		RXNORM_19AB_200406F	167	RXNORM	IN	167	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:09 PM	168		ENG	2020-04-06		RXNORM_19AB_200406F	168	RXNORM	IN	168	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:08 PM	173		ENG	2020-04-06		RXNORM_19AB_200406F	173	RXNORM	IN	173	Rxnorm	February 05, 2024		
02:03:47 PM	2020-04-27 09:04:13 PM	178		ENG	2020-04-06		RXNORM_19AB_200406F	178	RXNORM	IN	178	Rxnorm	February 05, 2024		

We can check all the data like wise .

- For the verification of row counts before and after the transformation we can check the log file which have logs for each file for the entire processing As shown below

```

94 1710148036256,"Processing file: RXNCONSO.RRF
95 "
96 1710148040120,"The count of RXNCONSO.RRF before transformation is: 1133065
97 "
98 1710148047748,"The RXNCONSO.RRF is transformed into CSV and stored into transformed_files in S3
99 "
100 1710148048407,"The count of RXNCONSO.RRF after transformation is: 1133065
101 "
102 1710148048407,"Processing file: RXNCUI.RRF
103 "
104 1710148048544,"The count of RXNCUI.RRF before transformation is: 30046
105 "
106 1710148048772,"The RXNCUI.RRF is transformed into CSV and stored into transformed_files in S3
107 "
108 1710148048799,"The count of RXNCUI.RRF after transformation is: 30046
109 "
110 1710148048799,"Processing file: RXNCUICHANGES.RRF
111 "
112 1710148048821,"The count of RXNCUICHANGES.RRF before transformation is: 153
113 "
114 1710148048862,"The RXNCUICHANGES.RRF is transformed into CSV and stored into transformed_files in S3
115 "
116 1710148048862,"The count of RXNCUICHANGES.RRF after transformation is: 153
117 "
118 1710148048862,"Processing file: RXNDOC.RRF
119 "
120 1710148048935,"The count of RXNDOC.RRF before transformation is: 3445
121 "
122 1710148049125,"The RXNDOC.RRF is transformed into CSV and stored into transformed_files in S3
123 "
124 1710148049126,"The count of RXNDOC.RRF after transformation is: 3445
125 "
126 1710148049126,"Processing file: RXNREL.RRF
127 "
128 1710148061108,"The count of RXNREL.RRF before transformation is: 7154306
129 "
130 1710148108526,"The RXNREL.RRF is transformed into CSV and stored into transformed files in S3

```