




Step 1: navigate to [Console Home | Console Home | us-east-1 \(amazon.com\)](#).




Step 2: Create a S3 Bucket



Step 3: Create a folder.



Step 4: Create a lambda function.



Step 5: Trigger the bucket

Lambda > Functions > rxnormdunction2

rxnormdunction2

Throttle Copy ARN Actions

✓ The trigger s3bucketexnorm was successfully added to function rxnormdunction2. The function is now receiving events from the trigger.

▼ Function overview Info

Export to Application Composer Download

Diagram Template

rxnormdunction 2

Layers (1)

S3

+ Add trigger

+ Add destination

Description

-

Last modified

1 minute ago

Function ARN

arn:aws:lambda:us-east-1:471112841584:function:rxnormdunction2

Function URL Info

-

Step 6: Deploy the code

✓ Successfully updated the function rxnormdunction2.

Code Test Monitor Configuration Aliases Versions

Code source Info

Upload from

File Edit Find View Go Tools Window Test Deploy

Go to Anything (Ctrl-P)

Environment

rxnormdunction2 - /

lambda_function.py

```
1 import json
2 import boto3
3 from zipfile import ZipFile
4 import pandas as pd
5 from io import BytesIO
6 from datetime import datetime
7 import os
8
9 s3_client = boto3.client('s3')
10
11 def lambda_handler(event, context):
12     try:
13         # Extract bucket and key information from the S3 event
14         bucket = event['Records'][0]['s3']['bucket']['name']
15         key = event['Records'][0]['s3']['object']['key']
16
17         if key.endswith('.xlsx'):
18             handle_excel_upload(bucket, key)
19         elif key.endswith('.zip'):
20             handle_zip_upload(bucket, key)
```

Step 7: Upload the files

Upload succeeded
View details below.

Upload: status

Close

The information below will no longer be available after you navigate away from this page.

Summary

Destination

s3://s3bucketexnorm

Succeeded

2 files, 238.6 MB (100.00%)

Failed

0 files, 0 B (0%)

Files and folders

Configuration

Files and folders (2 Total, 238.6 MB)

Find by name

< 1 >

Name	Folder	Type	Size	Status	Error
RxNorm_full...	-	application/...	238.5 MB	Succeeded	-
RxNorm_He...	-	application/...	24.7 KB	Succeeded	-

Step 8: Check the uploaded files.

Amazon S3 > Buckets > s3bucketexnorm > fileuploading/

fileuploading/

Copy S3 URI

Objects

Properties

Objects (9) Info

Copy S3 URI

Copy URL

Download

Open

Delete

Actions

Create folder

Upload

Find objects by prefix

Show versions

< 1 >

<input type="checkbox"/>	Name	Type	Last modified	Size	Storage class
<input type="checkbox"/>	RXNATOMARCHIVE.txt	txt	March 10, 2024, 23:59:24 (UTC-07:00)	58.9 MB	Standard
<input type="checkbox"/>	RXNCONSO.txt	txt	March 10, 2024, 23:59:33 (UTC-07:00)	138.5 MB	Standard
<input type="checkbox"/>	RXNCUI.txt	txt	March 10, 2024, 23:59:34 (UTC-07:00)	2.1 MB	Standard
<input type="checkbox"/>	RXNCUICHANGES.txt	txt	March 10, 2024, 23:59:34 (UTC-07:00)	17.9 KB	Standard
<input type="checkbox"/>	RXNDOC.txt	txt	March 10, 2024, 23:59:34 (UTC-07:00)	271.8 KB	Standard
<input type="checkbox"/>	RXNREL.txt	txt	March 11, 2024, 00:00:27 (UTC-07:00)	645.2 MB	Standard
<input type="checkbox"/>	RXNSAB.txt	txt	March 11, 2024, 00:00:33 (UTC-07:00)	10.5 KB	Standard
<input type="checkbox"/>	RXNSAT.txt	txt	March 11, 2024, 00:01:12 (UTC-07:00)	622.4 MB	Standard
<input type="checkbox"/>	RXNSTY.txt	txt	March 11, 2024, 00:01:20 (UTC-07:00)	26.1 MB	Standard

Step 9: Check the cloudwatch.

CloudWatch > Log groups > /aws/lambda/rxnordfunction2 > 2024/03/11/[\$LATEST]6da460e6006c440a9ca19395104e4cd5

Log events

Actions ▼

Start tailing

Create metric filter

You can use the filter bar below to search for and match terms, phrases, or values in your log events. [Learn more about filter patterns](#)

Q Filter events

1m1hLocal timezone ▼Display ▼

	Timestamp	Message
		No older events at this moment. Retry
	2024-03-10T23:58:07.547-07:00	INIT_START Runtime Version: python:3.12.v20 Runtime Version ARN: arn:aws:lam...
	2024-03-10T23:58:10.104-07:00	START RequestId: d6daf805-1385-402f-a95d-bf108ad59b05 Version: \$LATEST
	2024-03-10T23:58:10.105-07:00	.zip file uploaded: RxNorm_full_02052024.zip
	2024-03-10T23:58:12.653-07:00	Transformation Started
	2024-03-10T23:59:20.482-07:00	File Name Before Conversion After Conversion
	2024-03-10T23:59:20.482-07:00	0 RXNATOMARCHIVE 371368 371368
	2024-03-10T23:59:20.482-07:00	1 RXNCONSO 1133065 1133065
	2024-03-10T23:59:20.482-07:00	2 RXNCUI 30046 30046
	2024-03-10T23:59:20.482-07:00	3 RXNCUICHANGES 153 153
	2024-03-10T23:59:20.482-07:00	4 RXNDOC 3445 3445
	2024-03-10T23:59:20.482-07:00	5 RXNREL 7154306 7154306
	2024-03-10T23:59:20.482-07:00	6 RXNSAB 13 13
	2024-03-10T23:59:20.482-07:00	7 RXNSAT 7222404 7222404
	2024-03-10T23:59:20.482-07:00	8 RXNSTY 461874 461874
	2024-03-10T23:59:23.695-07:00	RXNATOMARCHIVE.txt is saved
	2024-03-10T23:59:33.593-07:00	RXNCONSO.txt is saved
	2024-03-10T23:59:33.730-07:00	RXNCUI.txt is saved
	2024-03-10T23:59:33.777-07:00	RXNCUICHANGES.txt is saved
	2024-03-10T23:59:33.887-07:00	RXNDOC.txt is saved
	2024-03-11T00:00:32.361-07:00	RXNREL.txt is saved
	2024-03-11T00:00:32.409-07:00	RXNSAB.txt is saved
	2024-03-11T00:01:18.312-07:00	RXNSAT.txt is saved
	2024-03-11T00:01:19.887-07:00	RXNSTY.txt is saved
	2024-03-11T00:01:19.887-07:00	Transformation Completed